

Viewing Subject Details Client 1.0.62

The search results from your query are presented both on the Summary tab in pie charts and on the Search Results tab in a table. Each record in the table corresponds to one subject (patient) and displays the collection(s), study(ies), and image series associated with that subject.

From the Search Results tab, you can view subject details, view all image series associated with the subject, and add image series to your cart for download. You can also view the DICOM header for a selected image series. In the future, you will be able to view the image.

To view subject details

1. **Search** for studies that match criteria that interest you.
The Summary tab refreshes to show your search results in pie chart form.
2. Click the **Search Results** tab.
The search results that match your criteria appear sorted by Subject ID.

Cart	Collection ID	Subject ID	Studies	Series
>	Head-Neck Cetuximab-Demo	0522c0001	4	13
>	Head-Neck Cetuximab-Demo	0522c0002	5	16
>	Head-Neck Cetuximab-Demo	0522c0003	5	16
>	Head-Neck Cetuximab-Demo	0522c0009	5	16
>	Head-Neck Cetuximab-Demo	0522c0012	3	11
>	Head-Neck Cetuximab-Demo	0522c0013	5	16
>	Head-Neck Cetuximab-Demo	0522c0014	7	14
>	Head-Neck Cetuximab-Demo	0522c0015	4	32
>	Head-Neck Cetuximab-Demo	0522c0017	4	26
>	Head-Neck Cetuximab-Demo	0522c0027	3	12

You can now **add** all studies associated with a subject to your cart at the same time or drill down to an individual study and add it to your cart instead.

3. Click a **Subject ID** link.
All studies associated with this subject appear. In this example, this subject is only associated with one study, and only one image series is associated with that study.

Collection ID	Subject ID	Studies	Series
Head-Neck Cetuximab-Demo	0522c0001	4	13
0522c0001	Neck*HeadNeckPETCT	3 Series	Study UID: ...54903251
	Aug 23, 1999*	4 Series	Study UID: ...2 28655
	Dec 07, 1999*	3 Series	Study UID: ...88515740
	Dec 07, 1999*	3 Series	Study UID: ...78294515

✔ The red asterisk indicates that the date, which is the date that this series was first available at TCIA, has been offset to protect private health information.

4. Click the arrow to the left of the cart icon to view details about the selected study.
All of the image series associated with that study appear.

Cart	Description	Modality	Manufacturer	Images	Series	Viewers	DICOM
>	SCOUT	CT	GE MEDICAL SYSTEMS	2	...01956153	Q	☰
>	C/A/P W/O CONTRAST	CT	GE MEDICAL SYSTEMS	123	...19889059	Q	☰
>	Recon 2: C/A/P W/O CONTRAST	CT	GE MEDICAL SYSTEMS	123	...33324740	Q	☰

You can now

- [Add an image series to your cart](#)
- View the series description, modality, manufacturer, number of images, and series ID
- [View thumbnails](#) of the image series
- [View an animation](#) of the image series
- [View DICOM data](#) for this image series